

# FlipMate® Filtration Systems



## Get quick filtration with ease

- Simply screw on filtration assembly, add an additional digestion cup, and apply vacuum
- Sample is ready in just seconds

**What's included:** FlipMate 50 systems include 200 filtration assemblies, 400 digestion cups (50 mL), and 200 caps. FlipMate 100 systems include 100 filtration assemblies, 200 digestion cups (100 mL), and 100 caps.



Membrane type	Pore size (µm)	Sample size	Complete FlipMate System		FlipMate Filtration Assembly Only	
			Cat. no.	Qty/pk	Cat. no.	Qty/pk
<b>FlipMate 50 systems</b>						
PES filter with prefilter	0.45	50 mL	SC0601C	200	SC0601	100
PES filter with prefilter, white caps	0.45	50 mL	SC0601CW	200	—	—
Glass fiber filter	1.5	50 mL	SC0602C	200	SC0602	100
Glass fiber filter, white caps	1.5	50 mL	SC0602CW	200	—	—
PES filter without prefilter	0.45	50 mL	SC0603C	200	SC0603	100
PES filter with PTFE prefilter†	0.45	50 mL	SC0604C	200	SC0604	100
Glass fiber filter with PTFE prefilter†	0.7	50 mL	SC0605C	200	SC0605	100
White glass fiber filter with PTFE prefilter†	0.7	50 mL	SC0605CW	200	—	—
Acid-washed glass fiber filter without prefilter	0.7	50 mL	—	—	SC0606	100
PVDF filter†	0.45	50 mL	—	—	SC0607	100
<b>FlipMate 100 systems</b>						
PES filter with prefilter	0.45	100 mL	SC0301C	100	SC0301	50
Glass fiber filter	1.5	100 mL	SC0302C	100	SC0302	50
PES filter with prefilter	0.2	100 mL	SC0303C	100	SC0303	50
Acid-washed glass fiber filter without prefilter	0.7	100 mL	SC0304C	100	SC0304	50
Acid-washed glass fiber filter without prefilter	>10	100 mL	SC0306C	100	SC0306	50
PES filter with PTFE prefilter†	0.45	100 mL	SC0308C	100	SC0308	50

†Certified for trace metals

Get full filtering in seconds—just 4 easy steps!

1. Bring sample to volume in appropriate digestion cup.
2. Screw FlipMate® assembly onto cup with sample and tighten securely.
3. Thread an additional digestion cup onto FlipMate assembly and tighten.
4. Invert system and apply vacuum.

## FlipMate® Vacuum Manifold

Filter up to 12 samples simultaneously

Description	Catalog number
PVC vacuum manifold, 12-place	SC0620



### Environmental Express

2345A Charleston Regional Parkway,  
Charleston, South Carolina 29492

Phone: 843-881-6560  
Toll Free: 800-343-5319  
Web: EnvExp.com  
E-mail: info@EnvExp.com

**EnvExp.com**