

GhostWipe® Wipes



Simplify your surface sample collection for lead and beryllium content

- Premoistened in deionized water and ready for sample collection
- Just wipe surface you want to test, digest, and analyze results
- Wipes breakdown completely and won't clog equipment
- Premoistened wipes comply with ASTM E 1792 for lead or beryllium in dust

Description	Catalog number	Qty/pk
GhostWipe premoistened wipes, 5.9" x 5.9" (15 x 15 cm), individually wrapped	SC4210	500
	SC4250	1000
GhostWipe dry wipes, 8" x 10" (20 x 25 cm)	SC4280	100

C1010 Dust sampling template, isolates a 3.9" x 3.9" (10 x 10 cm) area. 250/pk
C1012 Dust sampling template, isolates a 12" x 12" (30 x 30 cm) area. 100/pk
UC475-NL Ultimate Cup™ digestion cup, 50 mL, with natural linerless cap. 500/pk

Lead Wipe Standard in Cup

- Make certain of your accuracy
- GhostWipe® wipe spiked with in situ dust
- Wipes are loaded with known amounts of lead from 175 to 250 mg
- Each sample comes with sample ID# and certificate of analysis



Description	Catalog number	Qty/pk
Lead wipe standard in 50-mL cup	SC4252	10
Lead wipe standard in 100-mL cup	SC4256	10

Lead-Based Paint Standards

- Verify results with these standards
- Ground and sieved for homogeneity
- Each lot is digested and certified against NIST-traceable standards with both total and 3050B values reported

Description	Catalog number	Qty/pk
Lead-based paint chips, 0.1% Pb	HP-CRM-PC-A	40 grams



Place wipe in digestion cup and digest—wipes break down completely

Environmental Express

2345A Charleston Regional Parkway,
Charleston, South Carolina 29492

Phone: **843-881-6560**

Toll Free: **800-343-5319**

Web: **EnvExp.com**

E-mail: **info@EnvExp.com**

EnvExp.com