

Surface Sampling

GhostWipes from Environmental Express Simplify Sample Collection and Sample Prep for Lead or Beryllium in Surface Dust.


This sturdy wiping material moistened with DI water holds together even on the roughest wiping surfaces. In the lab, GhostWipes readily and completely dissolve during the digestion process providing more complete dispersion of analytes and more uniform recoveries. Complete breakdown of wipes eliminates fibrous residue that can clog expensive instruments. Each wipe is preserved in a sealed packet.

GhostWipe Options:

- Our original GhostWipe is **ASTM E 1792** compliant for Lead or Beryllium in Dust

Wipe size 5.9" x 5.9" (15cm x 15cm)	Catalog # SC4210	Pack of 500
	Catalog # SC4250	Pack of 1,000
- Our GhostWipe for Be is **ASTM D 7707** compliant for Beryllium in Dust

Wipe size 2" x 2.75" (5.1cm x 7cm)	Catalog # SC4210BE	Pack of 500
	Catalog # SC4250BE	Pack of 1,000



Take a sample with GhostWipe.



Digest as usual in acid solution.



GhostWipes totally dissolve during digestion for easy analysis.



Lead-based Paint Standard

These chips have been ground and sieved for homogeneity. Each lot is digested and certified for lead content using NIST traceable standards. Two digestion procedures are used and both total and 3050B values are reported.

- | | |
|-------------------------------------------|-----------------------|
| Lead-based Paint Chips, 0.1% Pb, 40 grams | Catalog # HP-CRM-PC-A |
| Lead-based Paint Chips, 0.5% Pb, 40 grams | Catalog # HP-CRM-PC-B |

Related Products:

- | | | |
|----------------------------------------------------|------------------|-------------|
| Dry Ghost Wipes in Individually Sealed Packets | Catalog # SC4280 | Case of 100 |
| Dust Sampling Template, 3.9" x 3.9" (10cm x 10cm) | Catalog # C1010 | Pack of 250 |
| Dust Sampling Template, 12" x 12" (30.5 c 30.5 cm) | Catalog # C1012 | Pack of 100 |
| Methamphetamine Sampling Template, 5cm (ID) | Catalog # C1014 | Pack of 200 |
| Sample Transport Vessel | Catalog # SC475 | Pack of 500 |

Lead Wipe Standard in Screw-Cap Rigid-walled Cup (SC475)

The lead wipe standard is a GhostWipe spiked with in-situ dust. Using an NIST traceable solid sample, known amounts of lead are placed on a wipe and each sample is individually labeled with the exact amount of lead the wipe contains. Each spike sample contains 175-250 micrograms of lead per wipe and come with a certificate of analysis showing the Sample ID# and the starting concentration of the dust.

Catalog # SC4252 Pack of 10

ENVIRONMENTAL EXPRESS

Call 800.745.8218 or 843.576.1147 ■ www.environmentalexpress.com

