

Organic Chemistry: Field Preservation Kits

Extend hold times on samples as much as 14 days with 5035 Field Sampling Kits from Environmental Express.

These convenient, economical, easy-to-use kits contain pre-labeled vials, reagents, certifications and instructions for soil sampling for volatile compounds by EPA Method 5035.

- All reagents and vials are certified with chromatograms
- No field weighing required
- Easy-to-follow instructions furnished with each kit
- Each kit is delivered in a convenient plastic rack that keeps vials together and upright

Volume discount pricing is available on all kits.

Custom kits are available. Contact our Technical Sales Department with your specifications.

All 5035 Field Preservation Kits Include:

- 5-gram coring device
- 40mL vial for dry weight sample
- 40mL weighed amber vial
- Convenient carrying kit
- Appropriate amount of Methanol in Snip and Pour Tube or Economical Vial
- Instructions for use

5035 Field Preservation Kits Sodium, Bisulfate/Methanol:

Two 40mL weighed vials containing 5mL NaHSO₄ with stir bars and Methanol tube or vial.

Description	Catalog #
Premium Kit, 5mL MeOH in PTFE tube	M3003
Economy Kit, 5mL MeOH in economical vial	M3013
Premium Kit, 10mL MeOH in PTFE tube	M3001
Economy Kit, 10mL MeOH in economical vial	M3011

5035 Field Preservation Kits, Organic-free Water and Methanol:

Two 40mL weighed vials with 5mL organic-free water with stir bars and Methanol tube or vial.

Description	Catalog #
Premium Kit, 5mL MeOH in PTFE tube	M3002
Economy Kit, 5mL MeOH in economical vial	M3012
Premium Kit, 10mL MeOH in PTFE tube	M3004
Economy Kit, 10mL MeOH economical vial	M3014

Methanol Snip and Pour PTFE tubes are sold in packs of 20.

Volume	Catalog #
5mL Methanol Snip	M3125
10mL Methanol Snip	M3127
15mL Methanol Snip	M3126
25mL Methanol Snip	M3128



Convenient Methanol Snips in Premium Field Kits eliminate the need to weigh Methanol in the field.

ENVIRONMENTAL EXPRESS

Call 800.343.5319 or 843.881.6560 • www.environmentalexpress.com

