

# StableWeigh™ TDS Vessels

## Eliminate cumbersome crucibles

- Save time—no more washing, drying, and weighing of hundreds of crucibles or dishes
- Preconditioned and preweighed—we've done the vessel prep work for you
- Faster cooling time—no more waiting for an hour
- Disposable—eliminate risk of carryover contamination
- Safer—stop safety hazards from broken crucibles

**Simplify your Total Dissolved Solids (TDS) analysis!** Made of proprietary materials, these revolutionary StableWeigh vessels will certainly replace the ceramic and glass crucibles you use in your lab today.

The name says it all—StableWeigh vessels are **weight stable**, even at the method-required temperature of 180°C. Each bag is predried and preweighed with the tare weight printed on the bag so it's **ready to use** right out of the box.

Think of the **time savings**—these disposable vessels remove all the preparatory drying and weighing steps prior to analysis, and the cleaning steps afterward.

Plus, you'll get a **better analysis**. Our innovative vessels weigh much less than an 80-gram crucible, so the weight ratio of vessel is closer to ending analysis result.

### Why Switch to StableWeigh™ Vessels?

- ✓ Preweighed with weight printed on the vessel
- ✓ Weight ratio of vessel is closer to ending analysis result
- ✓ Quicker to cool after heating
- ✓ Disposable—never have to wash
- ✓ Holds more volume than a typical dish

Description	Catalog number
StableWeigh TDS vessels; pack of 100	<b>TDS100</b>
Easy-use filling station for StableWeigh TDS vessels, 1 position	<b>TDS100F</b>
Easy-use filling station for StableWeigh TDS vessels, 6 position	<b>TDS600F</b>
Weighing stand for use on balance	<b>TDS200B</b>
Modular metal rack; holds multiple vessels for use in oven	<b>TDS200R</b>
Filtration kit: 100 TDS vessels and 100 washed and dried filters	<b>TDS111</b>
Starter kit: 100 TDS vessels, 100 washed and dried filters, weighing stand, oven rack, 6-position filling station, and static control system	<b>TDS222</b>

**Find out more about how StableWeigh™ vessels can improve your TDS analysis. Contact us today!**

#### Environmental Express

2345A Charleston Regional Parkway,  
Charleston, South Carolina 29492

Phone: **843-881-6560**

Toll Free: **800-343-5319**

Web: **EnvExp.com**

E-mail: **info@EnvExp.com**

 **EnvExp.com**